yac

Blaise Thompson

Introduction

Exercise

Using ya

yaq Instrument Shop Training

Blaise Thompson

University of Wisconsin-Madison

June 17, 2022




A Typical Instrument

Introduction

Evercise

Using yac





A Typical Instrument

Introduction

Exercis

Using ya

typical instrumental software...


- monolithic
 - hard to develop—need everything to work for anything to work
 - any one piece of the software can crash everything
 - hard to reuse pieces in other instruments
- inflexible
 - lots of "baked in" assumptions about the particular hardware attached
 - takes a long time to change experimental approach

... limits experimentalist freedom!



Evercies

Using ya





Exercis

Using ya

distributed

- each daemon can be developed separately from everything else
- no piece of software can crash everything
- possibility of multiclient or networked access

portable

- reusable daemons
- benefit from existing ecosystem rather than "reinventing" hardware support
- use just what you need
- multilingual (in theory)



Exercis

Using ya

Guiding principle: simplify client development as much as possible

YOU create clients



The yaq Approach

Blaise Thompson

Introduction

Using ya

self describing

- daemon tells the client about itself in a very structured way
- ▶ good documentation: https://yaq.fyi

traits

- enforce consistency between similar daemons where possible
- optional, extensible



Exercis

Using ya

simplify timing control

- any daemon can be asked for "busy" state
- architecture removes annoyance of blocking hardware interfaces



Exercise

Using ya

using your package manager, install yaqc

- client = yaqc.Client(host="192.168.1.1",
 port=38002)
- ▶ client.get_position
- ▶ client.set_position
- ▶ client.busy

Try comparing across all machines!



Exercise

Using ya

using your package manager, install yaqd-control

- yaqd scan --host 192.168.1.1 --start 38000
 --stop 38005
- ▶ yaqd status



Exercise

using your package manager, install yaqc-qtpy

graphical client based on traits


Test Hardware

Blaise Thompson

Introduction

Exercise

Using ya



fluorescein in ethanol fluorescence using Python, take an emission slice for a given excitation wavelength

- numpy linspace
- for loop
- matplotlib.pyplot



Using yaq

yaq doesn't have "built in" orchestration software

- write scripts
- create experiment-specific GUIs
- utilize existing projects focused on orchestration

the best choice depends on your instrument

- scripts are flexible and lightweight—rapid development, few users
- experiment-specific GUIs are great for mature experiments with many users
- can introduce as much or as little sophistication as needed

importantly you can choose a combination of many options!



Blaise Thompson

Introduction

Exercis

Using yaq



https://blueskyproject.io/

